Infrared Thermometer

User Manual

Model:

☐ A202C ☐ A202D ☐ A202E ☐ A202F

PLEASE NOTE:

THE MEDICAL DEVICE MUST BE USED ACCORDING TO USER MANUAL.

Shenzhen Aeon Technology Co., Ltd.

A301, Building A, Donghua Industrial Park, No. 5003, Bao'an Avenue, Sanwei Community,
Hangcheng Street, Bao'an District, 518126 Shenzhen, PEOPLE'S REPUBLIC OF CHINA.

Tel: +86-755-86182155

Customer Service E-mail: market@aeon-med.com

Website: www.aeon-med.com

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

E-mail: shholding@hotmail.com

File No.: A202D series IT-UM V 1.2 20260123

Content

1	Safety	4
1.1	Important safety instruction read before use	4
1.2	Cautions and Warnings	4
1.3	Definitions and Symbols	6
2	Introduction	7
2.1	Intended use.....	7
2.1.1	Intended user	7
2.1.2	Intended patient population	7
2.1.3	Medical condition.....	7
2.2	Contraindications.....	7
2.3	Intended patient population	8
2.4	Intended Users	8
2.5	Restrictions of Use	8
2.6	Display screen	8
2.7	Packing list	8
2.8	Product Features.....	9
3	Operation	9
3.1	Check.....	9
3.2	Install or replace batteries.....	9
3.3	Starting up	9
3.4	Placement	10
3.5	Mode switch	10
3.6	Memory mode and deletion.....	10
3.7	Unit switch	11
3.8	Sound switch settings	11
4	How to use your A202D series Infrared Thermometer.....	11
4.1	Body temperature measurement	11
4.2	Object temperature measurement	12
5	Errors and Troubleshooting	13
6	Care, Cleaning and Maintenance.....	13
7	Specification	14
8	Clinical accuracy	15
9	Disposal of device	15
10	Warranty	16
11	Manufacturer's EMC Declaration	16

Table 1: Device diagram of A202D Series infrared thermometer





Model	A202C	A202D
Annotated diagram		
Model	A202E	A202F
Annotated diagram		

Table 2: Description of the number in the product diagram above

No.	Description	No.	Description	No.	Description	No.	Description
1	Display	2	Power/Measurement Button	3	Forehead cover	4	Infrared Probe
5	Mode Switch	6	Label	7	Battery Cover	8	Disposable lens filter

1 Safety

1.1 Important safety instruction read before use

The following basic safety precautions should always be taken.

- (1) Close supervision is necessary when the thermometer is used by, on, or near children, handicapped persons or invalids.
- (2) Use the thermometer only for the intended use described in this manual.
- (3) Do not use the thermometer if it is not working properly, or if it has suffered any damage.

1.2 Cautions and Warnings

- (1) As with any thermometer, proper technique is crucial to obtaining accurate temperature readings. Please read this manual thoroughly and carefully before using.
- (2) Always operate the thermometer in an operating temperature range +5°C to +40°C, and relative humidity 15 to 95%.
- (3) Always store the thermometer in a cool and dry place -25°C to +55°C, and relative humidity 15% to 95%.
- (4) Although the thermometer works when battery indication appearing, the thermometer will not work.
- (5) Remove the batteries if stored for a long period of time.
- (6) The batteries should be kept away out of children's reach. If they are swallowed, promptly see a doctor for help.
- (7) The device requires no calibration. The product has been calibrated before the factory inspection.
- (8) When you find that the product does not measure correctly, please change the equipment to measure or contact your doctor and then contact the salesperson of that equipment for recalibration needs.
- (9) The device contains no user serviceable parts.
- (10) The user must check that the equipment functions safely and see that it is in proper working condition before being use.
- (11) The manufacturer does not require such preventive inspections by other persons.
- (12) No modification of this equipment is allowed. Do not disassemble the thermometer.
- (13) The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- (14) Manufacturer will provide circuit diagrams, component part lists, descriptions, calibration instructions to assist to service personnel in parts repair.
- (15) Do not clean or maintenance the device is in use. Do not touch the lens.
- (16) Store the thermometer in a cool and dry location. Avoid direct sunlight. Avoid dropping



the thermometer, if it happens and you think the thermometer may be damaged, please contact customer services immediately.

- (17) The effects of lint, dust, light (including sunlight), etc. when measuring.
- (18) Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.
- (19) This thermometer is not intended to substitute for a consultation with your physician.
- (20) This thermometer and the subject must remain in a stable environment for at least 30 minutes before measuring the temperature.
- (21) Please consult your doctor if the infrared thermometer shows elevated temperature (yellow or red backlights).
- (22) It takes at least 30 minutes for the infrared thermometer to warm from the minimum storage temperature between uses until it is ready for intended use.
- (23) The effects of degraded sensors and electrodes, or loosened electrodes, that can degrade performance or cause other problems.
- (24) Warnings regarding significant RISKS of reciprocal interference posed by ME EQUIPMENT during specific investigations or treatments.
- (25) Information on potential electromagnetic or other interference and advice on how to avoid or minimize such interference.
- (26) Warning statement that addresses the HAZARDS that can result from unauthorized modification of the ME EQUIPMENT.
- (27) Prolonged use Probe overheating, may burn the skin.
- (28) Do not use in an environment with high electromagnetic intensity, otherwise the equipment will be electrocuted and cause non-functioning, such as MRI.
- (29) Do not prevent the battery from being swallowed by kid in places where they can reach.
- (30) Do not use the device if the reading is incorrect, may be caused by low battery power.
- (31) Please remove the battery from the battery compartment to avoid overheating, when the device is not used for a long time.
- (32) Exceeding the expected service life (24 months) and inaccurate device measurements.
- (33) Please do not put the device near the fire or heat source to avoid the battery causing an explosion.
- (34) If the device is exposed to fluid ingress and is used before the cleaning fluid has completely dried, there is a risk of the device overheating potentially causing a burn to the user or the ear canal of the patient.
- (35) It cannot be used in patients with inflammation of the auditory canal or tympanic membrane.
- (36) It cannot be used in patients with obstruction of the external auditory canal.
- (37) It cannot be used in patients with hyperemia or earwax.

- (38) Please don't aim the source of infrared light at eyes.
- (39) Before measuring forehead temperature, make sure that the forehead is dry.
- (40) ASTM laboratory accuracy requirements in the display range of 37 to 39 °C (98 to 102 °F) for IR thermometers is $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$), whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$).
- (41) In order that purchasers may identify products conforming to requirements of this specification, producers and distributors may include a statement of compliance in conjunction with their name and address product labels or associated printed materials, or both, such as invoices, sales literature, and the like. The following statement is suggested: "This infrared thermometer meets requirements established in ASTM Standard (E 1965-98) Full responsibility for the conformance of this product to the standard is assumed by (name and address of producer or distributor)." In the event one or more provisions of this standard are not met, a cautionary statement shall be included.
- (42) Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.

1.3 Definitions and Symbols

Table 3: Description of symbol

Symbol	Description	Symbol	Description
	Type BF Equipment		Indicates a carry that contains unique device identifier information
	Indicates the manufacturer's serial number so that a specific medical device can be identified		Indicates the item is a medical device
	Information of manufacturer, including name and address		Indicates the date when the medical device was manufactured
	Indicates the temperature limits to which the medical device can be safely exposed		Indicates the range of humidity to which the medical device can be safely exposed
	Indicates the range of atmospheric pressure to which the medical device can be safely exposed		Follow the instructions for use

Symbol	Description	Symbol	Description
	CE mark, and the identification no. of notified body is 0123		Information of EU authorized representative
	For its disposal, this product must be sent to separate collection facilities for recovery and recycling		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences
<i>Note</i>	The important information you should know.	IP22	Anti-dust & Anti-water class

2 Introduction

2.1 Intended use

A202D Infrared Thermometer is intended for the intermittent measurement the human body temperature from forehead or ear canal measurement at home, clinics and hospital in people of all ages.

A202C and A202E is intended for the intermittent measurement the human body temperature from ear measurement at home, clinics and hospital in people of all ages. A202C must be measured with disposable lens filter of AEON.

A202F is intended for the intermittent measurement the human body temperature from forehead measurement at home, clinics and hospital in people of all ages.

2.1.1 Intended user

The Infrared Thermometer is intended for use by medical staff and lay persons.

2.1.2 Intended patient population

The Infrared Thermometer is intended to be used for people of any color and ages.

2.1.3 Medical condition

The Infrared Thermometer designed for measuring body temperature is intended to be used in clinical institutions and home environments.

2.2 Contraindications

Do not use for patients with wounds on their foreheads or ear canal.

Do not use for monitor patient continuously.

2.3 Intended patient population

The Infrared Thermometer is intended to be used for people of any color and ages.


2.4 Intended Users

The Infrared Thermometer is intended for use by medical staff and lay persons.

2.5 Restrictions of Use

The A202D Series infrared thermometer is clinically proven to produce accurate temperature measurements. However, please be advised that the accuracy cannot be ensured when the thermometer is not clean. Check that the probe is clean before taking a measurement.

2.6 Display screen

No.	Description	No.	Description
1	Memory symbol	2	Temperature display value
3	Calibration mode symbol	4	Sound switch symbol
5	Low battery symbol	6	Unit symbol
7	Decimal point	8	Ear mode symbol
9	Forehead mode symbol	10	0-36 months mode symbol
11	36 months mode symbol		

2.7 Packing list

Name of Articles	Quantity
Instruction manual, including warranty card	1
2 batteries, AAA 1.5V	1
Main engine	1
Disposable lens filter of 20 pcs	1

※Product packaging should contain the items described above. In case of any shortage, please contact Shenzhen Aeon Technology Co., Ltd. or the agent distributor in time.

2.8 Product Features

- (1) Simple and easy to use.
- (2) Automatically power off if left idle.
- (3) Memory function allows you to recall previous results up to 25 previous results.
- (4) Three colors of backlight, color visible of fever (red or yellow) and measurement in progress (green).

3 Operation

3.1 Check

Please check the packing case carefully before unpacking. In case of any damage found, please contact the carrier immediately. Open the package correctly, take out the infrared thermometer and other components from the case with care, and check them one by one against the packing list. When the equipment is moved a different environment, the difference in temperature or humidity may lead to condensation to it, in which case no use is allowed before condensation disappears.

3.2 Install or replace batteries

The first step after unpacking is to install the battery. The battery holder is on the back of the infrared thermometer. The battery installation steps are as follows:

- (1) Open the battery cover.
- (2) Insert the battery and keep +, - poles of the battery in line with +, - poles of the battery holder.
- (3) Close the battery cover.

Note:

Please use 2 AAA alkaline batteries.

Please do not mix-use old and new batteries if batteries are the same type.

Remove the batteries when the thermometer will not be used for a long time.

The Low battery symbol appearing on the screen indicates that the batteries run out; please replace the batteries before measuring.

3.3 Starting up

Press the Power button, then the backlight of the display comes on, the LCD is shown in full screen and the equipment enters the waiting state for measurement; at this time, the backlight goes out, put the probe on the forehead or ear canal, and the equipment gets ready for measurement. If there is no operation for about 10 seconds, the equipment will shut down

automatically.












3.4 Placement

Place the thermometer at a distance within 1-3cm from the center of the forehead or ear canal. If the eyebrow area is covered by hair, sweat or dirt, please clean it in advance to improve the accuracy of the reading. Keep the thermometer and forehead or ear still when measuring, as movement will negatively affect the temperature reading.

3.5 Mode switch

The device is an adjusted mode infrared thermometer. Toggle the slide switch to select the Forehead, ear mode (including 0-36 months mode and 36 months or more mode)/Calibration mode. Forehead or ear mode, i.e. operating mode, the measured result represents the equivalent temperature of the reference measurement part. The calibration mode, test mode, represents the temperature directly measured by the sensor and is used to verify the accuracy of the laboratory.

Table 4: Displayed mode of A202D Series infrared thermometer

Model	Calibration Mode	Forehead Mode		Ear Mode	
		0-36 Months	36 + Months	0-36 Months	36 + Months
A202D					
A202C, A202E		No such mode.			
A202F				No such mode.	

3.6 Memory mode and deletion

In shutdown mode, press and hold the power button, the screen displays "MEM" symbol, which flashes along with the symbol "MEM". Press the power button again to display the

number of memory groups 1 icon, and about 1 second later, show the display memory value, each time you press the power button, the memory value will be displayed. 25 groups will be stored in total, and when it is full.

In shutdown mode, remove the battery, and re-install it after 5 seconds, the memory value will be deleted.

3.7 Unit switch

In shutdown mode, press and hold the power button, it first displays "MEM", press and hold the button until "----°C" or "----°F" appears, press the power button to free convert between °C to °F, select °C or °F units and press and hold the power button to ensure that the units are saved correctly. It defaults to the current setting after rebooting.

3.8 Sound switch settings

In shutdown mode, press and hold the power button, it first displays "MEM", then "----°C" or "----°F" appears, last it displays "on + Sound switch symbol" or "off + Sound switch symbol", release the button at this time, it enters the sound setting mode; press the power button to change the sound switch once, press the power button to free convert between "on + Sound switch symbol" to "off + Sound switch symbol", a "bi" will be heard when the sound is enabled; press the power button again to change the sound, and there is no "bi" heard as the sound is disabled. Select on or off sound and press and hold the power button to ensure that the sound switch is saved correctly, It defaults to the current setting after rebooting.

4 How to use your A202D series Infrared Thermometer

4.1 Body temperature measurement

- (1) Toggle the “Mode switch” to the age precision mode  or  through each age group.

Note: you must select an age in order to take a measurement.

- (2) Choice your measurement side, forehead or ear canal. A202C and remove the forehead cover of A202D fit the probe snugly into the ear canal.

Note: Replace disposable hygienic lens filter with every use A202C to maintain hygiene and accuracy. Earwax and dirt on the tip of the lens filter, even in a microscopic layer, can create an obstacle between the sensor and the infrared heat emitted by the ear and impact the accuracy of the measurement. External factors like poor probe replacement, absent, defective, soiled disposable lens filter, used lens filter or dirty lens filter may influence ear temperature. Must use the company's

equipped with lens filter, so as to avoid the use of other models of lens filter caused may influence ear temperature.

- (3) A202C attach a new lens filter by pushing the thermometer probe straight into the lens filter inside the box and then pulling out.
- (4) Push and release the Power button. Align the thermometer probe to the position between eyebrows within 1-3cm from the forehead (suit for A202D forehead mode and A202F), or fit the probe snugly into the ear canal (suit for A202D ear mode, A202C and A202E); press the power button, when a "bi" is heard about 1s later (no "bi" heard if the sound is disabled), it means that the accurate body temperature measurement has been taken, with the result is shown on the LCD screen.

Age range	Green Normal temperature	Yellow Elevated temperature	Red High temperature
 0-36 Months	< 37.3°C (< 99.1°F)	≥ 37.3 — $\leq 37.9^{\circ}\text{C}$ (≥ 99.1 — $\leq 100.3^{\circ}\text{F}$)	> 37.9°C (> 100.3°F)
 36+ Months	< 37.5°C (< 99.5°F)	≥ 37.5 — $\leq 37.9^{\circ}\text{C}$ (≥ 99.5 — $\leq 100.3^{\circ}\text{F}$)	> 37.9°C (> 100.3°F)

- (5) Temperature displays.
- (6) The equipment will shut down automatically if there is no button press for (10±2) seconds.

Note:

- Since the measurement results are affected by sweat, oil and the surrounding environment, the measured temperature is for reference only.
- Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with the Infrared thermometer during or after breastfeeding, because the skin temperature may then be lower than the internal body temperature.
- It is advisable to take the temperature in the same ear, since the values measured in the left and right ear may be different.
- In the case of patients whose auditory canal is curved, the temperature values measured may deviate owing to anatomical influences.
- The ear should not be obstructed, for example by excessive amounts of ear wax.
- If the measured temperature is < 32.0°C (89.6°F) or > 43.0°C (109.4°F), the display will show with LO or HI symbol.

4.2 Object temperature measurement


- (1) Toggle the slide switch to the calibration mode, press the power button to turn on the thermometer, align the thermometer sensor to the object to be measured, at this time,

quickly adjust the appropriate distance within 1-3cm from the object surface, press the power button, when a "bi" is heard about 1S later (no "bi" heard if the sound is disabled), it means that the surface temperature has been measured, with the result displayed on the LCD screen.

- (2) Leave it idle for a few seconds the unit will switch off automatically.

5 Errors and Troubleshooting

Table 5: Troubleshooting

Error Message	Possible Reason	Resolution
No display On the LCD pane.	The battery has run out. Incorrect battery polarity.	Replace the battery. Please note: The (+) side of the battery must face upwards.
Measurement not possible (or an abnormal value is displayed).	The thermometer is not ready.	Wait until the °C symbol is displayed.
An abnormal temperature value is displayed.	The probe tip is dirty or damaged. Did you hear the beep after pressing the measurement button	Clean the probe tip or get it repaired. Wait until you hear the beep before removing the thermometer from the ear or forehead
	Temperature taken isn't within typical human temperature range outside the measuring Range (32.0°C - 43.0°C or 89.6°F - 109.4°F). Lo = temperature < 32.0°C (89.6°F). Hi = temperature > 43.0°C (109.4°F).	Check that the probe tip is clean and that the thermometer is properly placed on the forehead.
	The battery is low or too low to take correct temperature measurement.	Insert new batteries.
	Ambient temperature is not within the allowed operating temperature range (5°C - 40°C or 41°F - 104°F) or is changing too fast.	Allow the temperature to remain for 30 minutes in a room where the temperature is between 5°C - 40°C or 41°F - 104°F.

6 Care, Cleaning and Maintenance

The infrared thermometer has been calibrated before delivery from the factory. If you have any doubt about its accuracy, please contact the after-sales service. We recommend technical inspection for the measurement every two years, and it is required to comply with the

applicable national regulations of the local place. Technical inspection for the measurement may be carried out by the government agencies in charge or by authorized fee-paying maintenance services.

Never put the thermometer under a running tap or submerge it into water. Clean the thermometer shell and sensor with alcohol swab or cotton ball dipped in alcohol (95%), and be careful to avoid liquid entering the thermometer. Do not use abrasive cleaners.

7 Specification

Table 6: Specification of A202D Series infrared thermometer

Device Name	Infrared Thermometer
Dimensions(L×W×H) and Weight (including batteries)	A202C, A202E: (143×38×49) mm, about (78±2) g; A202D, A202F: (143×38×51.5) mm, about (85±2) g
Anti-electric Shock Type	Internally powered equipment
Anti-electric Shock Degree	Type BF equipment
EMC Type	Group I Class B
Enclosure Degree of ingress protection	IP22
Battery	2 × AAA, 1.5V alkaline batteries
Number of keys	2 keys, mode switch, measurement/power button
Temperature unit	Default °C or Default °F
Resolution	0.1°C or 0.1°F
Measurement range (Forehead or ear mode)	32.0°C - 43.0°C (89.6°F to 109.4°F)
Measurement accuracy (Forehead or ear mode)	± 0.2°C (± 0.4°F)
For forehead or ear mode T Indicator light range	1) 0-36 months Green: $T < 37.3^{\circ}\text{C}$ (99.1°F) Yellow: 37.3°C (99.1°F) $\leq T \leq 37.9^{\circ}\text{C}$ (100.3°F) Red: $T > 37.9^{\circ}\text{C}$ (100.3°F)
	2) 36+ months Green: $T < 37.5^{\circ}\text{C}$ (99.5°F); Yellow: 37.5°C (99.5°F) $\leq T \leq 37.9^{\circ}\text{C}$ (100.3°F); Red: $T > 37.9^{\circ}\text{C}$ (100.3°F)
Voice	Power on and ready for working: a short beep; Measurement finished:

	1 long beep equal or lower than 37.9°C (100.3°F), 6 short beeps higher than 37.9°C (100.3°F)
Memory	25 groups
Operating Environment	Temperature: +5°C to +40°C (+41°F to +104°F); Humidity: 15% to 95% RH (non-condensing); Air Pressure: 700hPa-1060hPa
Transport and Storage Environment	Temperature: -25°C to +55°C (-13°F to +131°F); Humidity: 15% to 95% RH (non-condensing); Air Pressure: 700hPa-1060hPa
Automatic power-off	(10±2) s
Software identify No. and Version	A202C: A202C Infrared thermometer embedded software, V1.1
	A202D: A202D Infrared thermometer embedded software, V1.1
	A202E: A202E Infrared thermometer embedded software, V1.1
	A202F: A202F Infrared thermometer embedded software, V1.1

8 Clinical accuracy

The below table shows statistic distribution of an invasive controlled desaturation study, which guided by cause 201.102 of ISO80601-2-56, Guideline for evaluating and documenting clinical Accuracy in human subjects.

Table 7: Clinical Accuracy

Model	Operating mode	Reference body site	Clinical Bias Δcb	Limits of Agreement LA	Clinical Repeatability σr
A202C	Ear	Armpit	-0.0363	0.4072	0.1423
A201D	Forehead	Armpit	-0.0906	0.5065	0.1663
	Ear	Armpit	-0.0594	0.4008	0.1526
A202E	Ear	Armpit	-0.0594	0.4008	0.1526
A202F	Forehead	Armpit	-0.0906	0.5065	0.1663

9 Disposal of device

Adhere to the applicable regulations when disposing of the device. This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove

the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

10 Warranty

The intended service life of the product is 5 years. The warranty period of the product is 1 year from the date of sale (see the package for manufacturing date).

Any damage caused by improper use, battery leakage, failure to operate in accordance with the requirements or transfer of the thermometer to other users is out of warranty. We will not provide free warranty service for the failure caused by the user's reasons as follows: Repair service beyond warranty will be charged accordingly.

- (1) Failure caused by unauthorized disassembly or refitting of the product.
 - (2) Failure caused by dropping during use or handling.
 - (3) Failure due to lack of proper maintenance.
 - (4) Failure to operate in accordance with the correct instructions in the operating manual.
- Failure caused by failing to operate correctly as instructed by the user manual, etc.

After-sale service unit: Shenzhen Aeon Technology Co., Ltd.

11 Manufacturer's EMC Declaration

- (1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
- (2) Guidance and manufacturer's declaration - Electromagnetic Emissions and Immunity.

Table 8 - Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions - CISPR 11	Group 1
RF emissions - CISPR 11	Class B
Harmonic emissions - IEC 61000-3-2	Not application
Voltage fluctuations/ flicker emissions - IEC 61000-3-3	Not application

Table 9 - Guidance and manufacturer's declaration - Electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 - Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Not applicable
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV.; line(s) to earth: ±2 kV.	Not applicable

	100 kHz repetition frequency	
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle; At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle; And 70% 25/30 cycles Single phase: at 0; 0% 300 cycles	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms, 6Vrms (in ISM and amateur radio bands); 80% Am at 1kHz	Not applicable
Radiated RF - IEC61000-4-3	10 V/m; 80 MHz – 2,7 GHz; 80 % AM at 1 kHz	10 V/m; 80 MHz – 2,7 GHz; 80 % AM at 1 kHz
Proximity magnetic fields IEC 61000-4-39	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 10 - Guidance and manufacturer's declaration - Electromagnetic Immunity							
	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity test Level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380–390	TETRA 400	Pulse Modulation 18 Hz	1,8	0.3	27
	450	430–470	GMRS 460, FRS 460	FM ±5kHz deviation 1 kHz sine	2	0.3	28
	710	704–787	LTE Band 13, 17	Pulse modulation 217 Hz	2	0.3	9
	745						
	780						
	810	800–960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band1,3,4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400–2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100–5800	WLAN 802.11a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5500						
	5785						

Table 11 - Guidance and manufacturer's declaration - electromagnetic Immunity		
Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz	CW	8
134,2 kHz	Pulse modulation a 2,1 kHz	65 b
13,56 MHz	Pulse modulation a 50 kHz	7,5 b
a) The carrier shall be modulated using a 50% duty cycle square wave signal.		
b) r.m.s., before modulation is applied.		

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Dział pomocy technicznej: Shenzhen Aeon Technology Co., Ltd.