Infrared Thermometer

User Manual

Model:

☐ A200 ☐ A201 ☐ A202

PLEASE NOTE:

THE MEDICAL DEVICE MUST BE USED ACCORDING TO USER MANUAL.

Shenzhen Aeon Technology Co., Ltd.

A301, Building A, Donghua Industrial Park, No. 5003, Bao'an Avenue, Sanwei Community,
Hangcheng Street, Bao'an District, 518126 Shenzhen, PEOPLE'S REPUBLIC OF CHINA.

Tel: +86-755-86182155

Customer Service E-mail: market@aeon-med.com

Website: www.aeon-med.com

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

E-mail: shholding@hotmail.com

File No.: A201 series IT-UM V 1.9 20260123

Content

1	Safety	2
1.1	Important safety instruction read before use	2
1.2	Cautions and Warnings	2
1.3	Definitions and Symbols	4
2	Introduction	5
2.1	Intended use.....	5
2.1.1	Intended user.....	5
2.1.2	Intended patient population.....	5
2.1.3	Medical condition	5
2.2	Contraindications.....	6
2.3	Restrictions of Use	6
2.4	How does it work.....	6
2.5	Product Features	6
3	Operation	6
3.1	Displayed mode	6
3.2	Select the Temperature Unit	7
3.3	Turn ON/OFF the LED Light.....	7
3.4	Replacing the Battery	8
3.5	Temperature measurement.....	8
3.5.1	For forehead mode	8
3.5.2	For ear mode (Suitable for A201).....	8
3.5.3	For Calibration mode	9
3.6	Memory	10
3.6.1	Recalling the memory	10
3.6.2	Clear the memory.....	10
4	Troubleshooting and Maintenance.....	10
4.1	Maintenance	10
4.2	Troubleshooting.....	11
5	Specification	11
6	Clinical accuracy	13
7	Disposal of device.....	13
8	Warranty	13
9	Manufacturer's EMC Declaration	14

Table 1: Device diagram of A201 Series infrared thermometer

Model	A200	A201
Annotated diagram		
Model	A202	
Annotated diagram		

Table 2: Screen display of A201 Series infrared thermometer

Model	A200	A201
Screen Display		
Model	A202	
Screen Display		

Table 3: Description of the number in the product diagram above

No.	Description	No.	Description	No.	Description
1	Infrared Probe	2	Labeling	3	Battery Cover
4	Display	5	Measurement Button	6	Power Button
7	LED Indicator Light	8	Mode Switch	9	Forehead cover
10	Memory Button	11	Calibration mode indication	12	Memory indication
13	Temperature data	14	Temperature unit	15	Battery indication
16	Forehead mode indication	17	Ear mode indication		

1 Safety

1.1 Important safety instruction read before use

The following basic safety precautions should always be taken.

- (1) Close supervision is necessary when the thermometer is used by, on, or near children, handicapped persons or invalids.
- (2) Use the thermometer only for the intended use described in this manual.
- (3) Do not use the thermometer if it is not working properly, or if it has suffered any damage.

1.2 Cautions and Warnings

- (1) As with any thermometer, proper technique is crucial to obtaining accurate temperature readings. Please read this manual thoroughly and carefully before using.
- (2) Always operate the thermometer in an operating temperature range +5°C to +40°C, and relative humidity 15% to 95%.
- (3) Always store the thermometer in a cool and dry place -25°C to +55°C, and relative humidity 15% to 95%.
- (4) The device requires no calibration. The product has been calibrated before the factory inspection.
- (5) When you find that the product does not measure correctly, please change the equipment to measure or contact your doctor and then contact the salesperson of that equipment for recalibration needs.
- (6) The device contains no user serviceable parts.
- (7) The user must check that the equipment functions safely and see that it is in proper working condition before being use.

- (8) The manufacturer does not require such preventive inspections by other persons.
- (9) No modification of this equipment is allowed. Do not disassemble the thermometer.
- (10) The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- (11) Manufacturer will provide circuit diagrams, component part lists, descriptions, calibration instructions to assist to service personnel in parts repair.
- (12) Do not clean or maintenance the device is in use. Do not touch the lens.
- (13) Store the thermometer in a cool and dry location. Avoid direct sunlight. Avoid dropping the thermometer, if it happens and you think the thermometer may be damaged, please contact customer services immediately.
- (14) The effects of lint, dust, light (including sunlight), etc. when measuring.
- (15) Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.
- (16) This thermometer is not intended to substitute for a consultation with your physician.
- (17) This thermometer and the subject must remain in a stable environment for at least 30 minutes before measuring the temperature.
- (18) Please consult your doctor if the infrared thermometer shows elevated temperature (yellow or red backlights).
- (19) It takes at least 30 minutes for the infrared thermometer to warm from the minimum storage temperature between uses until it is ready for intended use.
- (20) The effects of degraded sensors and electrodes, or loosened electrodes, that can degrade performance or cause other problems.
- (21) Warnings regarding significant RISKS of reciprocal interference posed by ME EQUIPMENT during specific investigations or treatments
- (22) Information on potential electromagnetic or other interference and advice on how to avoid or minimize such interference.
- (23) Warning statement that addresses the HAZARDS that can result from unauthorized modification of the ME EQUIPMENT.
- (24) The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.
- (25) Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.
- (26) Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- (27) Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or


decreased electromagnetic immunity of this equipment and result in improper operation.”

- (28) Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Infrared Thermometer (model name: A200, A201, A202), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- (29) A list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g., shielded cable, load impedance) or specifically (e.g., by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).
- (30) If any: the performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term “ESSENTIAL PERFORMANCE” need not be used).
- (31) The maximum temperature appeared on the ME equipment may exceed 41°C but less than 43°C for some area like display and battery compartment, in which the temperature is a safe skin contact temperature complying the Table 24 (Allowable maximum temperatures for skin contact with ME EQUIPMENT APPLIED PARTS) in IEC 60601-1:2005+AMD1:2012+AMD2:2020, the patient population shall be patients elder than 1 year of age, duration of contact < 10 min. Suitable for A202.
- (32) Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.

1.3 Definitions and Symbols

Table 4: Description of symbol

Symbol	Description	Symbol	Description
	Type BF Equipment		Indicates a carry that contains unique device identifier information
	Indicates the manufacturer's serial number so that a specific medical device can be identified		Indicates the item is a medical device

Symbol	Description	Symbol	Description
	Information of manufacturer, including name and address		Indicates the date when the medical device was manufactured
	Indicates the temperature limits to which the medical device can be safely exposed		Indicates the range of humidity to which the medical device can be safely exposed
	Indicates the range of atmospheric pressure to which the medical device can be safely exposed		Follow the instructions for use
	CE mark, and the identification no. of notified body is 0123		Information of EU authorized representative
	For its disposal, this product must be sent to separate collection facilities for recovery and recycling		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences
<i>Note</i>	The important information you should know.	IP22	Anti-dust & Anti-water class

2 Introduction

2.1 Intended use

The A201 Infrared Thermometer is intended for the intermittent measurement the human body temperature from forehead or ear measurement at home, clinics and hospital in people of all ages. A200 and A202 is intended for the intermittent measurement the human body temperature from forehead measurement at home, clinics and hospital in people of all ages.

2.1.1 Intended user

The Infrared Thermometer is intended for use by medical staff and lay persons.

2.1.2 Intended patient population

The Infrared Thermometer is intended to be used for people of any color and ages.

2.1.3 Medical condition

The Infrared Thermometer designed for measuring body temperature is intended to be used in clinical institutions and home environments.

2.2 Contraindications

Do not use for patients with wounds on their foreheads or ear. Do not use for monitor patient continuously.

2.3 Restrictions of Use

The A201 Series infrared thermometer is clinically proven to produce accurate temperature measurements.

However, please be advised that the accuracy cannot be ensured when the thermometer is not clean. Check that the probe is clean before taking a measurement.

2.4 How does it work

The device is an adjusted mode infrared thermometer. The A201 Series infrared thermometer measures the infrared heat generated by the surface of the skin over the vessel and its surrounding tissue. The thermometer then converts it into a temperature value.

Note: The A201 Series infrared thermometer does not emit any infrared energy.




2.5 Product Features








- (1) Automatically power off if left idle.
- (2) Memory function allows you to recall previous results up to 25 previous results.
- (3) Three colors of backlight, color visible of fever (red or yellow) and measurement in progress (green).

3 Operation

3.1 Displayed mode

Table 5: Displayed mode of A201 series infrared thermometer

Mode	A200	A201	A202	Mark
Forehead Mode				The mode is used to measure the forehead temperature.

Mode	A200	A201	A202	Mark
Calibration Mode				The mode is used to measure the object temperature.
Ear Mode	No such mode.		No such mode.	The mode is used to measure the ear temperature.
Memory Mode				View Measured Values.

3.2 Select the Temperature Unit

The A201 Series infrared thermometer provides two measurement units used for indicating the body/object temperature, °C or °F, for your preferred selection.

- (1) In the off state, press and hold the measurement button with your finger until ‘---°C/°F’ displayed, then remove your finger.
- (2) Press the measurement button again, the unit will change to ‘---°F/°C’, each time you press the measurement button, the unit will change once, select the unit you need to show on the screen.
- (3) Then press the power button or let the thermometer turn off automatically (suitable for A200 and A201).
- (4) Then press and hold the power button thermometer until it turns off (suitable for A202).

3.3 Turn ON/OFF the LED Light



Only suitable for A200 and A202, not A201. The LED Light to help users placing the thermometer at the correct position.

A200: In the off state, press and hold the measurement button with your finger until ‘---°C/°F’ displayed, then remove your finger. Press the power button again, the display will change to ‘LED/OFF’, each time you press the measurement button, the display will change once, select the you need to show on the screen. Then press the power button or let the thermometer turn off automatically.

A202: In the off state, press and hold the measurement button with your finger until ‘LED’ displayed, then remove your finger. Press the measurement button again, the LED light change. Each time you press the measurement button, but the LED light changes once, according to your need for the LED light on/off, press and hold the measurement button again until the backlight goes off.

3.4 Replacing the Battery

The A201 series infrared thermometer comes with two 1.5V AAA alkaline batteries. The meter will display battery indication to alert you when the meter power is getting low, please follow the steps below to replace new batteries immediately.

	Remove the battery cover as the arrow direction accordingly.
	Remove the old batteries and replace with two new 1.5V AAA size alkaline batteries. Taking care to match the Positive (+) and Negative (-) indications.
	Close the battery cover as the arrow direction accordingly.

Note:

- (1) Although the thermometer works when battery indication appearing, the infrared thermometer will not work.
- (2) Remove the batteries if stored for a long period of time.
- (3) The batteries should be kept away out of children's reach. If they are swallowed, promptly see a doctor for help.

3.5 Temperature measurement

3.5.1 For forehead mode

- (1) Turn on the infrared thermometer by first pressing the power button. Press the mode switch to select forehead mode. Forehead mode indication is displayed on the screen.
- (2) Place the probe close to the forehead and take a measurement, making sure the probe is flat and close to the forehead, not at an angle, within 3 cm distance. Press the measurement button once, while there may be a short beeping sound.
- (3) The measurement takes approximately 1-2s, a long beep will sound when the measurement is finished and the temperature measurement value is displayed along with the backlight color. To take the next measurement, wait for an audible reminder at the end of the previous measurement.
- (4) Press the power button to turn off the unit, or leave it idle for a few seconds the unit will switch off automatically.

3.5.2 For ear mode (Suitable for A201)

- (1) Press the power button, turn on the unit and a beep is heard, remove the forehead cover of infrared thermometer. Ear mode indication is displayed on the screen. The infrared thermometer is now in ear mode. (Cover the forehead cover to enter the forehead mode, and the forehead mode indication is displayed on the screen).
- (2) If the person whose temperature is being taken is a child under one year old, make the ear canal straight by gently pulling the ear backward. Other ages of the person whose temperature was measured, make the ear canal straight by gently pulling the ear backward and upward.
- (3) Insert the probe into the auditory canal. For the temperature to be measured accurately in the ear canal, the tip of the probe must be oriented vertically towards the eardrum.
- (4) Press the Measurement button and wait until you hear a short beep and a long beep. As soon as you hear the long beep, remove the probe from the auditory canal. Read the temperature displayed on the screen.

Note:

- Since the measurement results are affected by sweat, oil and the surrounding environment, the measured temperature is for reference only.
- Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with the Infrared thermometer during or after breastfeeding, because the skin temperature may then be lower than the internal body temperature.
- It is advisable to take the temperature in the same ear, since the values measured in the left and right ear may be different.
- In the case of patients whose auditory canal is curved, the temperature values measured may deviate owing to anatomical influences.
- The ear should not be obstructed, for example by excessive amounts of ear wax.
- If the measured temperature is $< 32.0^{\circ}\text{C}$ (89.6°F) or $> 43.0^{\circ}\text{C}$ (109.4°F), the display will show with LO or HI symbol.

3.5.3 For Calibration mode

- (1) Switching the thermometer to calibration mode. A200 and A202: just need to move the mode switch to the upper end, at this time, the thermometer display calibration mode indication. A201: Press the power button and measurement button at the same time, wait until the calibration mode indication is displayed and then release.
- (2) This mode can measure the temperature of the object. Move the probe to within 3 cm of the surface of the object, and when ready, press the measurement button and take the measurement.
- (3) Read the result. The measurement result will be done in 1 second and one long beep

informs about the temperature measurement, after a while another short beep, confirming the saving of the result to the memory and readiness for the next measurement.

- (4) Press the power button to turn off the unit, or leave it idle for a few seconds the unit will switch off automatically.

3.6 Memory

3.6.1 Recalling the memory

The thermometer can store up to 25 recent temperature measurements.

- (1) In the off state, A200 and A201 is to press the measurement button once the display appears 'MEM', A202 is to press and hold the measurement button with your finger until the "MEM" disappears, then remove your finger.
- (2) When you press the measurement button once, the first times and test data will alternate on the display, the records that have been tested will be shown in reverse order (the latest result is shown first), and green, yellow or red background light will appear according to the memory reading.
- (3) When the number of measurement results exceeds 25, the oldest result cannot be found by the memory function because a new result will be added. When the last record is shown on the display, press the measurement button again to return to the latest record of the measurement result.
- (4) Exit Memory mode. Press the Power Button to turn off the device, or let the device sit idle for a few seconds and the device will automatically turn off.

3.6.2 Clear the memory

A200 and A201: In the off state, press the measurement button once the display appears 'MEM' THEN remove your finger. press and hold the power button until 'dEL' is shown on the display and release your finger, press and hold the power button until 'CLr' is shown on the display and release your finger. All memory data will be deleted.

A202: In the off state, remove the battery cover, then press and hold the measurement button with your hand, at the same time, remove a battery with the other hand, all the memory data will be deleted at this time.

4 Troubleshooting and Maintenance

4.1 Maintenance

Never put the thermometer under a running tap or submerge it into water. Clean the thermometer shell and sensor with alcohol swab or cotton ball dipped in alcohol (95%), and be careful to avoid liquid entering the thermometer. Do not use abrasive cleaners.

4.2 Troubleshooting

Table 6: Troubleshooting

Problem	Possible Reason	Resolution
No display On the LCD pane.	The battery has run out. Incorrect battery polarity.	Replace the battery. Please note: The (+) side of the battery must face upwards.
Measurement not possible (or an abnormal value is displayed).	The thermometer is not ready.	Wait until the °C symbol is displayed.
An abnormal temperature value is displayed.	The probe tip is dirty or damaged. Did you hear the beep after pressing the measurement button?	Clean the probe tip or get it repaired. Wait until you hear the beep before removing the thermometer from the ear or forehead
LO or HI symbol is displayed	The temperature measured is outside the measuring Range. LO-temperature < 32.0°C (89.6°F). HI-temperature > 43.0°C (109.4°F).	Check that the probe tip is clean and that the thermometer is properly placed on the forehead.
 Symbol is displayed	The battery has run out.	Replace the battery.
 Symbol is displayed.	The ambient temperature is outside the operating temperature range or is changing too fast.	To ensure accurate measurement, let the thermometer rest at operating temperature for 30 minutes prior to use.

5 Specification

Table 7: Specification of A201 Series infrared thermometer

Device Name	Infrared Thermometer
Dimensions(L×W×H) and Weight (including batteries)	A200: (170×44×54) mm, about (98±2) g; A201: (159×44×50) mm, about (82±2) g; A202: (143×38×40) mm, about (80±2) g

Anti-electric Shock Type	Internally powered equipment
Anti-electric Shock Degree	Type BF equipment
EMC Type	Group I Class B
Enclosure Degree of ingress protection	IP22
Battery	2×AAA, 1.5V alkaline batteries
Number of keys	A200: 3 keys, switch button, measurement button, power button; A201: 2 keys, measurement button, memory button; A202: 2 keys, switch button, measurement/power button
Temperature unit	Default °C or Default °F
Resolution	0.1°C or 0.1°F
Measurement range (Forehead or ear mode)	32.0°C- 43.0°C (89.6°F - 109.4°F)
Measurement accuracy (Forehead or ear mode)	±0.2°C (±0.4°F)
For forehead or ear mode T Indicator light range	Green: $T < 37.5^{\circ}\text{C}$ (99.5°F) Yellow: 37.5°C (99.5°F) $\leq T \leq 37.9^{\circ}\text{C}$ (100.3°F) Red: $T > 37.9^{\circ}\text{C}$ (100.3°F)
Voice	Power on and ready for working: a short beep; Measurement finished: 1 long beep equal or lower than 37.9°C (100.3°F), 6 short beeps higher than 37.9°C (100.3°F)
Memory	25 groups
Operating Environment	Temperature: +5°C to +40°C (+41°F to +104°F); Humidity: 15% to 95% RH (non-condensing); Air Pressure: 700hPa-1060hPa
Transport and Storage Environment	Temperature: -25°C to +55°C (-13°F to +131°F); Humidity: 15% to 95% RH (non-condensing); Air Pressure: 700hPa-1060hPa
Automatic power-off	A200: (30±2) s; A201, A202: (12±2) s
Software identify No.	A200: A200 Infrared thermometer embedded software, V1.2
	A201: A201 Infrared thermometer embedded software, V1.2
	A202: A202 Infrared thermometer embedded software, V1.2

6 Clinical accuracy

The below table shows statistic distribution of an invasive controlled desaturation study, which guided by cause 201.102 of ISO80601-2-56, Guideline for evaluating and documenting clinical Accuracy in human subjects.

Table 8: Clinical Accuracy

Model	Operating mode	Reference body site	Clinical Bias Δcb	Limits of Agreement LA	Clinical Repeatability σr
A200	Forehead	Axillary	-0.1163	0.3678	0.1057
A201	Forehead	Axillary	-0.1113	0.3859	0.0899
	Ear	Axillary	-0.0388	0.3911	0.0909
A202	Forehead	Axillary	-0.1250	0.3209	0.1315

7 Disposal of device

Adhere to the applicable regulations when disposing of the device. This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

8 Warranty

The intended service life of the product is 5 years. The warranty period of the product is 1 year from the date of sale (see the package for manufacturing date).

Any damage caused by improper use, battery leakage, failure to operate in accordance with the requirements or transfer of the thermometer to other users is out of warranty. We will not provide free warranty service for the failure caused by the user's reasons as follows: Repair service beyond warranty will be charged accordingly.

- (1) Failure caused by unauthorized disassembly or refitting of the product.
- (2) Failure caused by dropping during use or handling.
- (3) Failure due to lack of proper maintenance.
- (4) Failure to operate in accordance with the correct instructions in the operating manual.

Failure caused by failing to operate correctly as instructed by the user manual, etc.

After-sale service unit: Shenzhen Aeon Technology Co., Ltd.

9 Manufacturer's EMC Declaration

- (1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
- (2) Guidance and manufacturer's declaration - Electromagnetic Emissions and Immunity.

Table 9 - Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions - CISPR 11	Group 1
RF emissions - CISPR 11	Class B
Harmonic emissions - IEC 61000-3-2	Not application
Voltage fluctuations/ flicker emissions - IEC 61000-3-3	Not application

Table 10 - Guidance and manufacturer's declaration - Electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 - Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Not applicable
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV; line(s) to earth: ±2 kV. 100 kHz repetition frequency	Not applicable
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle; At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle; And 70% 25/30 cycles Single phase: at 0; 0% 300 cycles	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms, 6Vrms (in ISM and amateur radio bands); 80% Am at 1kHz	Not applicable
Radiated RF - IEC61000-4-3	10 V/m; 80 MHz – 2,7 GHz; 80 % AM at 1 kHz	10 V/m; 80 MHz – 2,7 GHz; 80 % AM at 1 kHz
Proximity magnetic fields IEC 61000-4-39	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m
NOTE U _T is the a.c. mains voltage prior to application of the test level.		

Table 11 - Guidance and manufacturer's declaration - Electromagnetic Immunity							
Radiated RF IEC61000-4-3	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity test Level (V/m)

(Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380–390	TETRA 400	Pulse Modulation 18 Hz	1,8	0.3	27
	450	430–470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704–787	LTE Band 13, 17	Pulse modulation 217 Hz	2	0.3	9
	745						
	780						
	810	800–960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band1,3,4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400–2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100–5800	WLAN 802.11a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5500						
	5785						

Table 12 - Guidance and manufacturer's declaration - electromagnetic Immunity

Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz	CW	8
134,2 kHz	Pulse modulation ^a 2,1 kHz	65 ^b
13,56 MHz	Pulse modulation ^a 50 kHz	7,5 ^b
a) The carrier shall be modulated using a 50% duty cycle square wave signal.		
b) r.m.s., before modulation is applied.		

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Modell	A200	A201
Kommentiertes Diagramm	Diagram of the A200 infrared thermometer. The front view shows the lens (1), display (2), and battery cover (3). The side view shows the trigger (4), power button (5), and other controls (6, 8).	Diagram of the A201 infrared thermometer. The front view shows the lens (1), display (2), and battery cover (3). The side view shows the trigger (4), power button (5/6), and other controls (8, 10). A separate view shows the lens cap (9).
Modell	A202	
Kommentiertes Diagramm	Diagram of the A202 infrared thermometer. The front view shows the lens (1), display (2), and battery cover (3). The side view shows the trigger (4), power button (5/6), and other controls (8).	

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Modell	A200	A201
Bildschirmanzeige	Diagram of the A200 LCD display. It shows icons for home (11), person (12), and temperature (14). The display shows four digits (13) and a battery level indicator (15). The text 'MEM' (12) is also visible.	Diagram of the A201 LCD display. It shows icons for home (11), person (12), and temperature (14). The display shows four digits (13) and a battery level indicator (15). The text 'MEM' (12) is also visible.
Modell	A202	
Bildschirmanzeige	Diagram of the A202 LCD display. It shows icons for home (11), person (12), and temperature (14). The display shows four digits (13) and a battery level indicator (15). The text 'MEM' (12) is also visible.	

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Dział pomocy technicznej: Shenzhen Aeon Technology Co., Ltd.